AimStrip® TANDEM MEASURING SYSTEM

LIPID PROFILE & GLUCOSE MEASURING SYSTEM

AimStrip[®] Tandem Meter is used for testing Cholesterol, Triglycerides and Blood Glucose levels in capillary whole blood.¹





Use with the AimStrip[®] brand test strips and control solutions - designed for exclusive use with the AimStrip[®] Tandem Meter.

LOOK AT THESE GREAT FEATURES!

Rapid

- Simply insert test, apply specimen & read results
- No specimen preparation required

POC Testing

- Handheld meter for POC testing
- ⇒ Easy to read LCD results

Minimal Training Time

User friendly meter makes testing simple

Zero Maintenance

- No servicing of meter is required
- → Meter has a 3-year warranty

Convenient

- Automatic code recognition by RFID-tag
- ⇒ Battery operated

CAP / CLSI Procedures

- → Already written, available as MS Word & PDF
- → Just add your information

CPT CODE: Total Cholesterol	.82465 QW*
CPT CODE: Triglycerides	.84478 QW*
CPT CODE: Blood Glucose	.82947 QW*

THE TESTS YOU WANT. THE ACCURACY YOU NEED!

DISTRIBUTED BY:



¹ See Package Insert for complete procedure. *CPT Code is assigned by Medicare. Check with each state for exact reimbursement.

Size (mm): 61 W x 109 D x 23 H

CLIA

WAIVED

Weight: 77.5g (including battery) ±1g



AimStrip[®] Tandem Measuring System

- AimStrip[®] Tandem Meter
- → Warranty Card

- ⇒ Log Book
- Store at 2° to 30°C







REF#

77301

AimStrip

12:03

RFID

(((...)))

CHOL

3

мем

TG

MULTIPLE

PATIENT USE

ТС

AimStrip® TANDEM MEASURING SYSTEM

CLIA WAIVED

REF# 77310

77410

77550

REF# 77312

77412

77512

AimStrip® Tandem TEST STRIPS + CONTROL SOLUTIONS For use with the AimStrip[®] Tandem Measuring System Only

\frown		
Ainst	AimStrip [®] Tandem	
Amotrip tandem	TEST STRIPS AimStrip [®] Tandem	SIZE
Comp R	TOTAL CHOLESTEROL	10 Strips
	🗢 10 Test Strips	
Ainstrip Tandam	→ 10 Capillary Rods	
	AimStrip [®] Tandem TRIGLYCERIDES	10 Strips
Canada	🗢 10 Test Strips	
	🗢 10 Capillary Rods	
Ainting Lander and the	AimStrip [®] Tandem	
50	BLOOD GLUCOSE	50 Strips
	Packed in 2 Bottles	
	🗢 25 Test Strips per Bottle	
STORE AT 8° TO 30°C	AimStrip [®] Tandem	
8 10 30 C	CONTROL SOLUTIONS	
		CIZE
	AimStrip [®] Tandem TOTAL CHOLESTEROL	SIZE 2 x 3mL
		2 X SITIL
	 ➡ Control Solution Level 1 (3mL) ➡ Control Solution Level 2 (3mL) 	
	AimStrip [®] Tandem	
	TRIGLYCERIDES	2 x 3mL
	Control Solution Level 2 (3mL)	
	AimStrip [®] Tandem	
	GLUCOSE	3 x 3.5mL
	➡ Control Solution Level 1 (3.5mL)	
	➡ Control Solution Level 2 (3.5mL)	AimStrip Tanders
6		And the second s
		Ginning

-DATE 12:03 • Ċ MEM RFID (((--))) тс

PRODUCT SPECIFICATIONS

Test	Measuring Principle	Spectroscopy	
Lile	Testing Types	Total Cholesterol, Triglycerides	
Pro	Testing Method	Enzymatic-colorimetric method	
Lipid Profile Test	Testing Range	Total Cholesterol: 100-400 mg/dL Triglycerides: 70-600 mg/dL	
	Testing Time	Approx 2 min	
	Calibration	Plasma - Equivalent	
	Sample Volume	5µL	
	HCT Range	30 - 55%	
Glucose Test	Measuring Principle	Electrochemical	
ose	Testing Types	Glucose	
iluce	Testing Range	20-600 mg/dL	
0	Testing Time	3 sec	
	Calibration	Plasma - Equivalent	
	Sample Volume	0.3µL	
	HCT Range	30 - 55%	
RFID	Standard Supported	ISO/IEC 15693	
	Nominal Read/Write Range	< 2cm	
	Antenna	Integrated	
	Frequency	13.56 MHz	
	Power Supply	(AAA) 1.5V x 2	
	Battery Lifetime	Approx 1,000 tests	
	PC Port	For cable connection between Meter and PC or thermal printer	

THE TESTS YOU WANT. THE ACCURACY YOU NEED!



www.germainelabs.com